

Anti-ABCC12 antibody

Catalog Number: 175833

Product name

Anti-ABCC12 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to ABCC12

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human mrp9. aa range:691-740

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000


Validations

Fig1.; Western Blot analysis of various cells using ABCC12 Polyclonal Antibody diluted at 1: 500