

SULT4A1 antibody



Catalog Number: 115833

Product name

SULT4A1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SULT4A1 Rabbit Polyclonal antibody. Positive WB detected in rat brain tissue, human brain tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 33kd

Preparation

This antibody was obtained by immunization of SULT4A1 recombinant protein (Accession Number: NM_014351). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IP, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

Validations



rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115833(SULT4A1 Antibody) at dilution of



1:600

IP Result of anti-SULT4A1 (IP:Catalog No:115833, 4ug; Detection:Catalog No:115833 1:300) with mouse brain tissue lysate 4000ug.