

SULT4A1 antibody



Catalog Number: 115833

Product name

ELISA, IP, WB

SULT4A1 antibody

Dilutions

Specificity

Recommended Dilution:

Human, Mouse, Rat; other species not tested.


WB: 1:200-1:2000

Antibody description

IP: 1:200-1:1000

SULT4A1 Rabbit Polyclonal antibody. Positive WB detected in rat brain tissue, human brain tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 33kd

Validations



Preparation

This antibody was obtained by immunization of SULT4A1 recombinant protein (Accession Number: NM_014351). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:115833(SULT4A1 Antibody) at dilution of

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality


Polyclonal

Ig Type

1:600

Rabbit IgG

Applications



IP Result of anti-SULT4A1 (IP:Catalog No:115833, 4ug; Detection:Catalog No:115833 1:300) with mouse brain tissue lysate 4000ug.