

Neurotrophin 3 antibody



Catalog Number: 113133

Product name

Neurotrophin 3 antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

Specificity


Human, Mouse, Rat; other species not tested.

IF: 1:20-1:200

Validations

Antibody description

Neurotrophin 3 Rabbit Polyclonal antibody.
Positive IP detected in SH-SY5Y cells. Positive WB detected in SH-SY5Y cells, human placenta tissue. Positive IHC detected in human brain tissue. Positive IF detected in SH-SY5Y cells. Observed molecular weight by Western-blot: 35 kDa




Preparation

This antibody was obtained by immunization of Neurotrophin 3 recombinant protein (Accession Number: NM_002527). Purification method: Antigen affinity purified.

SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:113133(NTF3 antibody) at dilution of 1:1000

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Immunohistochemical of paraffin-embedded human brain using Catalog No:113133(NTF3 antibody) at dilution of 1:100 (under 10x lens)


Applications

ELISA, WB, IHC, IF, IP

Dilutions

Recommended Dilution:


WB: 1:1000-1:10000



Neurotrophin 3 antibody


Catalog Number: 113133

Immunohistochemical of paraffin-embedded human brain using Catalog No:113133(NTF3 antibody) at dilution of 1:100 (under 40x lens)



IP Result of anti-NTF3 (IP:Catalog No:113133, 4ug; Detection:Catalog No:113133 1:800) with

SH-SY5Y cells lysate 1600ug.



Immunofluorescent analysis of (-20oc Ethanol) fixed SH-SY5Y cells using Catalog No:113133(NTF3 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)