

PPM1B antibody

Catalog Number: 114133

Product name

PPM1B antibody

IHC: 1:20-1:200

IF: 1:50-1:500


Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PPM1B Rabbit Polyclonal antibody. Positive IF detected in MCF-7 cells. Positive IHC detected in human breast cancer tissue, human skeletal muscle tissue. Positive WB detected in MCF7 cells. Observed molecular weight by Western-blot: 50-53 kDa

Validations




Preparation

This antibody was obtained by immunization of PPM1B recombinant protein (Accession Number: NM_177969). Purification method: Antigen affinity purified.

MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:114133(PPM1B antibody) at dilution of 1:400

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:114133(PPM1B Antibody) at dilution of 1:200 (under 10x lens)

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC, IF

Dilutions

Recommended Dilution:

WB: 1:200-1:2000




Immunohistochemistry of paraffin-embedded

PPM1B antibody

Catalog Number: 114133

human breast cancer tissue slide using Catalog No:114133(PPM1B Antibody) at dilution of 1:200 (under 40x lens)



Immunofluorescent analysis of (-20oc Ethanol) fixed MCF-7 cells using Catalog No:114133(PPM1B Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)