

# TIMM44 antibody



Catalog Number: 116133

## Product name

TIMM44 antibody

## Immunogen

[Human TIMM44 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TIMM44 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, A549 cells, Jurkat cells. Positive IF detected in A549 cells. Observed molecular weight by Western-blot: 45 kDa

## Preparation

This antibody was obtained by immunization of TIMM44 recombinant protein (Accession Number: NM\_006351). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116133(TIMM44 antibody) at dilution of 1:500



Immunofluorescent analysis of A549 cells, using TIMM44 antibody Catalog No:116133 at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).