

Anti-PKR Mouse antibody



Catalog Number: 169133

Product name

Anti-PKR Mouse antibody

Specificity

Human

Antibody description

PKR Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human EIF2AK2 (AA: 329-551) expressed in E. Coli.

Formulation

Purified antibody in PBS with 0.05% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Mouse IgG1


Applications

WB, FC

Validations




Western blot analysis using EIF2AK2 mAb against



human EIF2AK2 (AA)

Western blot analysis using EIF2AK2 mouse mAb against A431 (1), MCF-7 (2), PC-12 (3) cell lysate.



Flow cytometric analysis of A431 cells using EIF2AK2 mouse mAb (green) and negative control

Anti-PKR Mouse antibody



Catalog Number: 169133

Western blot analysis using EIF2AK2 mAb against HEK293 (1) and EIF2AK2 (AA

