

# Anti-PKR Mouse antibody



Catalog Number: 169133

## Product name

Anti-PKR Mouse antibody

## Specificity

Human

## Antibody description

PKR Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant fragment of human EIF2AK2 (AA: 329-551) expressed in E. Coli.

## Formulation

Purified antibody in PBS with 0.05% sodium azide.

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal


## Ig Type

Mouse IgG1


## Applications

WB, FC

## Validations




Western blot analysis using EIF2AK2 mAb against



human EIF2AK2 (AA)

Western blot analysis using EIF2AK2 mouse mAb against A431 (1), MCF-7 (2), PC-12 (3) cell lysate.



Flow cytometric analysis of A431 cells using EIF2AK2 mouse mAb (green) and negative control

# Anti-PKR Mouse antibody



Catalog Number: 169133

---

Western blot analysis using EIF2AK2 mAb against  
HEK293 (1) and EIF2AK2 (AA

