

# HOXB4 antibody



Catalog Number: 111533

## Product name

HOXB4 antibody

## Immunogen

[Human HOXB4 Recombinant protein \(GST tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

HOXB4 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 30-35 kDa

## Preparation

This antibody was obtained by immunization of HOXB4 recombinant protein (Accession Number: NM\_024015). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:111533(HOXB4 Antibody) at dilution of 1:600