

p120 Catenin antibody



Catalog Number: 113533

Product name

p120 Catenin antibody

WB: 1:200-1:2000

IP: 1:200-1:2000

Specificity

Human, Mouse, Rat; other species not tested.


IHC: 1:50-1:500

IF: 1:10-1:100

Antibody description

p120 Catenin Rabbit Polyclonal antibody. Positive WB detected in RAW264.7 cells, human brain tissue, mouse brain tissue, NIH/3T3 cells. Positive IP detected in mouse brain tissue. Positive IF detected in HepG2 cells. Positive IHC detected in human colon cancer tissue. Observed molecular weight by Western-blot: 105 kDa, 65-75 kDa

Validations



Preparation

This antibody was obtained by immunization of p120 Catenin recombinant protein (Accession Number: NM_001085468). Purification method: Antigen affinity purified.

RAW264.7 cells were subjected to SDS PAGE followed by western blot with Catalog No:113533(CTNND1 antibody) at dilution of 1:500

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications

ELISA, WB, IF, IP, IHC

Dilutions


Recommended Dilution:



Immunofluorescent analysis of HepG2 cells, using CTNND1 antibody Catalog No:113533 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit


p120 Catenin antibody

Catalog Number: 113533




IgG (red).

IP Result of anti-CTNND1 (IP:Catalog No:113533, 4ug; Detection:Catalog No:113533 1:500) with mouse brain tissue lysate 2640ug.



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using Catalog No:113533(CTNND1 Antibody) at dilution of 1:200 (under 10x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using Catalog No:113533(CTNND1 Antibody) at dilution of 1:200 (under 40x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).