

PRPF3 Polyclonal Antibody



Catalog Number: 164533

Product name

PRPF3 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to PRPF3

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human PRPF3 (NP_004689.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF, IP

Dilutions


WB 1:500 - 1:2000

IHC 1:50 - 1:200

IF 1:50 - 1:200


IP 1:50 - 1:100

Validations




Western blot - PRPF3 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using PRPF3 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunohistochemistry - PRPF3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using PRPF3 Antibody at dilution




of 1:100 (40x lens).

PRPF3 Polyclonal Antibody

Catalog Number: 164533


Immunohistochemistry - PRPF3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human gastric cancer using PRPF3 Antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - PRPF3 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human breast using PRPF3 Antibody at dilution of



1:100 (40x lens).


Immunohistochemistry - PRPF3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human liver cancer using PRPF3 Antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - PRPF3 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded mouse brain using PRPF3 Antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - PRPF3 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using PRPF3 Antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - PRPF3 Polyclonal Antibody

PRPF3 Polyclonal Antibody



Catalog Number: 164533

Immunohistochemistry of paraffin-embedded rat brain using PRPF3 Antibody at dilution of 1:100 (40x lens).