

# TEAD3 antibody

Catalog Number: 115933

---

## Product name

TEAD3 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

TEAD3 Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, mouse lung tissue. Observed molecular weight by Western-blot: 53kd

## Preparation

This antibody was obtained by immunization of TEAD3 recombinant protein (Accession Number: BC027877). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:115933(TEAD3 antibody) at dilution of 1:500