

# AKAP4 antibody

Catalog Number: 107933

## Product name

AKAP4 antibody

## Specificity

Human; other species not tested.

## Antibody description

AKAP4 Rabbit Polyclonal antibody. Positive WB detected in human testis tissue. Positive IHC detected in human breast cancer tissue. Observed molecular weight by Western-blot: 82 kDa

## Preparation

This antibody was obtained by immunization of AKAP4 recombinant protein (Accession Number: NM\_003886). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:107933(AKAP4 Antibody) at dilution of 1:1500



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:107933(AKAP4 Antibody) at dilution of 1:50 (under 40x lens)