

DVL2 Polyclonal Antibody



Catalog Number: 162233

Product name

DVL2 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to DVL2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human DVL2 (NP_004413.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


IHC, IF

Dilutions

IHC 1:50 - 1:200

IF 1:50 - 1:200

Validations



Immunofluorescence - DVL2 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using DVL2 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.