

NR1I3 antibody



Catalog Number: 113233

Product name

NR1I3 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

NR1I3 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HepG2 cells, L02 cells. Observed molecular weight by Western-blot: 45kd

Preparation

This antibody was obtained by immunization of NR1I3 recombinant protein (Accession Number: BC069651). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113233(NR1I3 antibody) at dilution of 1:300