

Human TMX2/TXNDC14/CGI-31 (His tag) recombinant protein



Catalog Number: 518233

General Information

Gene Name Synonym

Thioredoxin-related transmembrane protein 2,
Cell proliferation-inducing gene 26 protein,
Thioredoxin domain-containing protein 14, TMX2

Protein Construction

A DNA sequence encoding the Human
TMX2/TXNDC14/CGI-31 Met125-Lys296 is
expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL
method.

Stability

Store for up to 12 months at -20°C to -80°C as
lyophilized powder. Aliquots of reconstituted
samples are stable at $< -20^{\circ}\text{C}$ for 3 months.

Formulation

Lyophilized from a $0.2\ \mu\text{m}$ filtered solution of
20mM Tris, 150mM NaCl, pH8.0.

Usage Guide


Storage

Lyophilized protein should be stored at $< -20^{\circ}\text{C}$,
though stable at room temperature for 3 weeks.
Reconstituted protein solution can be stored at
 $4-7^{\circ}\text{C}$ for 2-7 days.

Reconstitution

Always centrifuge tubes before opening. Do not
mix by vortex or pipetting. It is not recommended
to reconstitute to a concentration less than
 $100\ \mu\text{g}/\text{ml}$. Dissolve the lyophilized protein in
ddH₂O. Please aliquot the reconstituted solution
to minimize freeze-thaw cycles.

SDS-PAGE



Human TMX2/TXNDC14/CGI-31 (His tag)
recombinant protein