

# Human TMX2/TXNDC14/CGI-31 (His tag) recombinant protein



Catalog Number: 518233

## General Information

### Gene Name Synonym

Thioredoxin-related transmembrane protein 2,  
Cell proliferation-inducing gene 26 protein,  
Thioredoxin domain-containing protein 14, TMX2

### Protein Construction

A DNA sequence encoding the Human  
TMX2/TXNDC14/CGI-31 Met125-Lys296 is  
expressed with a 6His tag at the N-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL  
method.

### Stability

Store for up to 12 months at -20°C to -80°C as  
lyophilized powder. Aliquots of reconstituted  
samples are stable at < -20°C for 3 months.

## Formulation

Lyophilized from a 0.2 µm filtered solution of  
20mM Tris, 150mM NaCl, pH8.0.

## Usage Guide


### Storage

Lyophilized protein should be stored at <-20°C,  
though stable at room temperature for 3 weeks.  
Reconstituted protein solution can be stored at  
4-7°C for 2-7 days.

### Reconstitution

Always centrifuge tubes before opening. Do not  
mix by vortex or pipetting. It is not recommended  
to reconstitute to a concentration less than  
100µg/ml. Dissolve the lyophilized protein in  
ddH<sub>2</sub>O. Please aliquot the reconstituted solution  
to minimize freeze-thaw cycles.

## SDS-PAGE



Human TMX2/TXNDC14/CGI-31 (His tag)  
recombinant protein