

# CHMP2A antibody



Catalog Number: 109233

## Product name

CHMP2A antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CHMP2A Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HepG2 cells. Observed molecular weight by Western-blot: 32 kDa

## Preparation

This antibody was obtained by immunization of CHMP2A recombinant protein (Accession Number: NM\_014453). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:109233(CHMP2A antibody) at dilution of 1:300