

# FGF17 antibody



Catalog Number: 110633

## Product name

FGF17 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

FGF17 Rabbit Polyclonal antibody. Positive IF detected in NIH/3T3 cells. Positive WB detected in MCF-7 cells, PC-3 cells. Observed molecular weight by Western-blot: 33 kDa

## Preparation

This antibody was obtained by immunization of FGF17 recombinant protein (Accession Number: NM\_003867). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IF, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

## Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:110633(FGF17 Antibody) at dilution of 1:600



Immunofluorescent analysis of NIH/3T3 cells using Catalog No:110633(FGF17 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)