

MEF2C antibody



Catalog Number: 112633

Product name

MEF2C antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MEF2C Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, mouse small intestine tissue, mouse testis tissue, rat lymph tissue, SH-SY5Y cells. Positive IHC detected in human heart tissue. Observed molecular weight by Western-blot: 65-68kd

Preparation

This antibody was obtained by immunization of MEF2C recombinant protein (Accession Number: BC026341). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:112633(MEF2C antibody) at dilution of 1:500



Immunohistochemistry of paraffin-embedded human heart slide using Catalog No:112633(MEF2C Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded human heart slide using Catalog No:112633(MEF2C Antibody) at dilution of 1:50