

# ST8SIA1 antibody



Catalog Number: 115633

## Product name

ST8SIA1 antibody

## Specificity

Human; other species not tested.

## Antibody description

ST8SIA1 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells, MDA-MB-453s cells. Observed molecular weight by Western-blot: 43 kDa

## Preparation

This antibody was obtained by immunization of ST8SIA1 recombinant protein (Accession Number: NM\_003034). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:115633(ST8SIA1 Antibody) at dilution of 1:1000