

DEDD Polyclonal Antibody



Catalog Number: 166633

Product name

DEDD Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to DEDD

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-318 of human DEDD (NP_001034801.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:100

Validations



Immunohistochemistry - DEDD Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat lung using DEDD antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - DEDD Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human esophagus using DEDD antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - DEDD Polyclonal Antibody

DEDD Polyclonal Antibody



Catalog Number: 166633

Immunohistochemistry of paraffin-embedded

mouse kidney using DEDD antibody at dilution of 1:100 (40x lens).