

TPI1 antibody



Catalog Number: 107633

Product name

TPI1 antibody

Specificity

Human; other species not tested.

Antibody description

TPI1 Mouse Monoclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells. Observed molecular weight by Western-blot: 27kd

Preparation

This antibody was obtained by immunization of TPI1 recombinant protein (Accession Number: NM_000365). Purification method: Caprylic Acid/Ammonium Sulfate Precipitation.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:107633(TPI1 antibody) at dilution of 1:500