

Anti-PFKFB3 Rabbit antibody



Catalog Number: 168633

Product name

Anti-PFKFB3 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

PFKFB3 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human PFKFB3

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC, IP

Dilutions


WB: 1/1000-1/5000

ICC/IF: 1/20

FC: 1/20

IP: 1/20

Validations



Western blot detection of PFKFB3 in K562 cell lysates using PFKFB3 Rabbit mAb(1:1000 diluted). Predicted band size:60kDa. Observed band size:60kDa.