

CX3CL1 antibody



Catalog Number: 109633

Product name

CX3CL1 antibody

Specificity

Human; other species not tested.

Antibody description

CX3CL1 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, human brain tissue. Positive IHC detected in human small intestine tissue. Observed molecular weight by Western-blot: 90-100 kDa

Preparation

This antibody was obtained by immunization of CX3CL1 recombinant protein (Accession Number: NM_002996). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:109633(CX3CL1 antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human small intestine slide using Catalog No:109633(CX3CL1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human small intestine slide using Catalog No:109633(CX3CL1 Antibody) at dilution of 1:50

CX3CL1 antibody



Catalog Number: 109633

(under 40x lens)