

# Mouse Angptl4/Farp/Fiaf/Ng27 (Fc tag) recombinant protein



Catalog Number: 519633

## General Information

### Gene Name Synonym

Angiopoietin-related protein 4,425O18-1,Angiopoietin-like protein 4,Fasting-induced adipose factor,Hepatic fibrinogen/angiopoietin-related protein,HFARP,Secreted protein Bk89,Angptl4,Farp, Fiaf, Ng27

### Protein Construction

A DNA sequence encoding the Mouse Angptl4/Farp/Fiaf/Ng27 Lys167-Ser410 is expressed with a Fc tag at the C-terminus.

### Organism

Mouse

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

### Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of PBS,pH7.4.

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

### SDS-PAGE



Mouse Angptl4/Farp/Fiaf/Ng27 (Fc tag)  
recombinant protein