

LRRC41 antibody



Catalog Number: 112333

Product name

LRRC41 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

LRRC41 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Positive IHC detected in human testis tissue, human liver tissue. Observed molecular weight by Western-blot: 80-90 kDa

Preparation

This antibody was obtained by immunization of LRRC41 recombinant protein (Accession Number: BC004953). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112333(LRRC41 antibody) at dilution of 1:100



Immunohistochemical of paraffin-embedded human testis using Catalog No:112333(LRRC41 antibody) at dilution of 1:50 (under 10x lens)