

Anti-FOLR1 / Folate Receptor 1 antibody



Catalog Number: 105333

Product name

Anti-FOLR1 / Folate Receptor 1 antibody

Immunogen

[Rat FOLR1 / Folate Receptor 1 \(His Tag\) recombinant protein](#)

Specificity

Rat FOLR1 / Folate Receptor 1

Antibody description

Rabbit polyclonal to FOLR1 / Folate Receptor 1

Preparation

Produced in rabbits immunized with purified, recombinant Rat FOLR1 / Folate Receptor 1 (rR FOLR1 / Folate Receptor 1; NP_598211.1; Met1-Met231). FOLR1 / Folate Receptor 1 specific IgG was purified by Rat FOLR1 / Folate Receptor 1 affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC-P

Dilutions

WB: 1-5 µg/ml


IHC-P: 0.1-2 µg/mL

Validations



Items	Lanes	A
Sample (whole cell lysate)		K562
Sample Volume (µg/lane)		30
Gel		13% SDS-PAGE reducing gel
Recommended Concentration		10 µg/ml
Secondary Antibody		Dylight 800 labeled Antibody to Rabbit IgG (H+L), at 1:5000 dilution.
		Developed using Odyssey imaging system.
Explanation		Predicted band size : 30 kDa Observed band size : 37 kDa

FOLR1 / Folate Receptor 1 Antibody, Rabbit PAb, Antigen Affinity Purified, Western blot



FOLR1 / Folate Receptor 1 Antibody, Rabbit PAb, Antigen Affinity Purified, Immunohistochemistry

Immunochemical staining of rat FOLR1 in rat kidney with rabbit polyclonal antibody (1 µg/mL, formalin-fixed paraffin embedded sections).