

transgelin/SM22 antibody



Catalog Number: 116333

Product name

transgelin/SM22 antibody

Immunogen

[Human TAGLN Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

transgelin/SM22 Rabbit Polyclonal antibody.
Positive WB detected in mouse colon tissue,
human colon tissue. Positive IHC detected in
human heart tissue, human colon cancer tissue.
Observed molecular weight by Western-blot:
19-22kd

Preparation

This antibody was obtained by immunization of
transgelin/SM22 recombinant protein (Accession
Number: NM_003186). Purification method:
Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH
7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



mouse colon tissue were subjected to SDS PAGE
followed by western blot with Catalog
No:116333(TAGLN Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded
human heart tissue slide using Catalog
No:116333(TAGLN Antibody) at dilution of 1:50
(under 10x lens)



transgelin/SM22 antibody



Catalog Number: 116333

Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog

No:116333(TAGLN Antibody) at dilution of 1:50 (under 40x lens)