

Human LIPF (His tag) recombinant protein



Catalog Number: 519333

General Information

Gene Name Synonym

Gastric Triacylglycerol Lipase, GL, Gastric Lipase, LIPF

Protein Construction

A DNA sequence encoding the Human LIPF Leu20-Lys398 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 25mM TrisHCl, 100mM glycine, 10%glycerol,pH 7.3.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human LIPF (His tag) recombinant protein