

FMR1NB antibody



Catalog Number: 110733

Product name

FMR1NB antibody

Specificity

human; other species not tested.

Antibody description

FMR1NB Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, A549 cells, SKOV-3 cells. Observed molecular weight by Western-blot: 66kd

Preparation

This antibody was obtained by immunization of FMR1NB recombinant protein (Accession Number: NM_152578). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:110733(FMR1NB antibody) at dilution of 1:1000