

ACTL6A Polyclonal Antibody



Catalog Number: 163733

Product name

ACTL6A Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to ACTL6A

Preparation

Antigen: A synthetic peptide of human ACTL6A

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000

IF 1:20 - 1:50

Validations



Western blot - ACTL6A Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using ACTL6A antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.