

C4orf7 antibody



Catalog Number: 108733

Product name

C4orf7 antibody

Immunogen

[Human C4orf7 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

C4orf7 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells. Observed molecular weight by Western-blot: 8 kDa

Preparation

This antibody was obtained by immunization of C4orf7 recombinant protein (Accession Number: NM_152997). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:1000

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:108733(C4orf7 Antibody) at dilution of 1:300