

# CTRP4 antibody

Catalog Number: 109733

## Product name

CTRP4 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CTRP4 Rabbit Polyclonal antibody. Positive WB detected in BxPC-3 cells. Positive IHC detected in human brain tissue. Observed molecular weight by Western-blot: 45kd

## Preparation

This antibody was obtained by immunization of CTRP4 recombinant protein (Accession Number: XM\_017017166). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



BxPC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:109733(C1QTNF4 antibody) at dilution of 1:1500



Immunohistochemical of paraffin-embedded human brain using Catalog No:109733(C1QTNF4 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human brain using Catalog No:109733(C1QTNF4 antibody) at dilution of 1:100 (under 40x lens)