

Human PLA2G1B/PLA2/PLA2A/PPLA2 (His tag) recombinant protein



Catalog Number: 519733

General Information

Gene Name Synonym

PLA2G1B/Phospholipase A2/Group IB
phospholipase A2/PLA2/PLA2A/PPLA2

Protein Construction

A DNA sequence encoding the Human
PLA2G1B/PLA2/PLA2A/PPLA2 Ala23-Ser148 is
expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL
method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM
TrisHCl, 150mM NaCl, 10% Glycerol, pH8.0.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human PLA2G1B/PLA2/PLA2A/PPLA2 (His tag)
recombinant protein