

Human FAM20C / DMP4 (His Tag) recombinant protein



Catalog Number: 501043

General Information

Gene Name Synonym

Dentin matrix protein 4; Golgi-enriched fraction casein kinase; Protein FAM20C

Protein Construction

A DNA sequence encoding the human FAM20C (NP_064608.2) (Met1-Arg584) was expressed with a polyhistidine tag at the C-terminus.

Organism

Human

Expression Host

Baculovirus-Insect Cells

QC Testing

Purity

> 85 % as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Leu 23

Molecular Mass

The recombinant human FAM20C consists 573 amino acids and predicts a molecular mass of 65.1 kDa.

Formulation

Lyophilized from sterile 20 mM Tris, 500 mM NaCl, pH 8.0, 10 % glycerol.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human FAM20C / DMP4 Protein (His Tag)