

NCOA5 antibody

Catalog Number: 113043

Product name

NCOA5 antibody

IP: 1:200-1:2000

IF: 1:10-1:100

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NCOA5 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells. Positive IP detected in HEK-293 cells. Positive FC detected in HEK-293 cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 66 kDa

Preparation

This antibody was obtained by immunization of NCOA5 recombinant protein (Accession Number: BC056872). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IF, IP, FC

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:113043(NCOA5 antibody) at dilution of 1:600



Immunofluorescent analysis of HeLa cells, using NCOA5 antibody Catalog No:113043 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit




IP Result of anti-NCOA5 (IP:Catalog No:113043,

NCOA5 antibody



Catalog Number: 113043

4ug; Detection:Catalog No:113043 1:800) with HEK-293 cells lysate 3000ug.



1×10^6 HEK-293 cells were stained with 0.2ug NCOA5 antibody (Catalog No:113043, red) and control antibody (blue). Fixed with 4% PFA blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1500.