

Bcl9 Polyclonal Antibody

Catalog Number: 165043

Product name

Bcl9 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to Bcl9

Preparation

Antigen: Recombinant protein of mouse Bcl9

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations




Western blot - Bcl9 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using Bcl9 antibody at 1:400 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.



Immunohistochemistry - Bcl9 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat brain using BCL9 antibody at dilution of 1:100



(40x lens).

Bcl9 Polyclonal Antibody



Catalog Number: 165043

Immunohistochemistry - Bcl9 Polyclonal Antibody

mouse brain using BCL9 antibody at dilution of 1:100 (40x lens).

Immunohistochemistry of paraffin-embedded