

# Anti-phospho-Rb (Thr252) Rabbit antibody



Catalog Number: 168043

## Product name

Anti-phospho-Rb (Thr252) Rabbit antibody

## Ig Type

Rabbit IgG

## Specificity

Human

## Applications

WB, IHC-P

## Antibody description

phospho-Rb (Thr252) Rabbit monoclonal antibody

## Dilutions

WB: 1/1000

## Preparation

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Thr252 of human Rb

IHC: 1/200

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA


## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Validations



Western blot detection of Rb (phospho-Thr252) in HL-60,U2OS,C6 using Rb (phospho-Thr252) Rabbit mAb(1:1000 diluted)