

# SLC17A7 Polyclonal Antibody



Catalog Number: 163443

---

## Product name

SLC17A7 Polyclonal Antibody

## Specificity

Mouse, Rat

## Antibody description

Polyclonal antibody to SLC17A7

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 461-560 of human SLC17A7 (NP\_064705.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:200

## Validations



Immunofluorescence - SLC17A7 Polyclonal Antibody

Immunofluorescence analysis of mouse brain using SLC17A7 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.