

# NUPL1 antibody



Catalog Number: 113443

---

## Product name

NUPL1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

NUPL1 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, HeLa cells, L02 cells, mouse brain tissue. Observed molecular weight by Western-blot: 65 kDa

## Preparation

This antibody was obtained by immunization of NUPL1 recombinant protein (Accession Number: NM\_014089). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:113443(NUPL1 antibody) at dilution of 1:500