

Anti-FGF11 antibody



Catalog Number: 176443

Product name

Anti-FGF11 antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to FGF11

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the internal region of human fgf11. aa range:21-70

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of NIH-3T3, HeLa cells using FGF-11 Polyclonal Antibody.. Secondary antibody (catalog#: HA1001) was diluted at 1:20000