

Anti-Lipocalin-2/NGAL (4C12) Mouse antibody



Catalog Number: 168443

Product name

Anti-Lipocalin-2/NGAL (4C12) Mouse antibody

Specificity

Human

Antibody description

Lipocalin-2/NGAL (4C12) Mouse monoclonal antibody

Preparation

Antigen: Recombinant protein of human NGAL.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Mouse IgG2b

Applications

ELISA

Validations



Indirect ELISA assay for anti-NGAL mouse mAb. Antigen coating concentration: 4ug/ml.