

# TAL1 antibody

Catalog Number: 115843

## Product name

TAL1 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

TAL1 Rabbit Polyclonal antibody. Positive IP detected in Raji cells. Positive WB detected in Raji cells, Jurkat cells. Observed molecular weight by Western-blot: 40-44 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_003189). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

## Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:115843(TAL1 antibody) at dilution of 1:1000



IP Result of anti-TAL1 (IP:Catalog No:115843, 4ug; Detection:Catalog No:115843 1:1000) with Raji cells lysate 2000ug.