

AHCYL1 antibody



Catalog Number: 107843

Product name

AHCYL1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

AHCYL1 Rabbit Polyclonal antibody. Positive IF detected in Hela cells. Positive WB detected in Jurkat cells, HeLa cells, mouse thymus tissue. Observed molecular weight by Western-blot: 60kd

Preparation

This antibody was obtained by immunization of AHCYL1 recombinant protein (Accession Number: NM_006621). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:107843(AHCYL1 antibody) at dilution of 1:1000



Immunofluorescent analysis of Hela cells, using AHCYL1 antibody Catalog No:107843 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).