

Mapatumumab reference antibody



Catalog Number: 177143

Product name

Mapatumumab reference antibody

Specificity

Human TNFRSF10A / TRAILR1 / CD261

Antibody description

Anti-TNFRSF10A / TRAILR1 / CD261 Antibody (Mapatumumab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG2SA

Applications


ELISA

Validations



SDS-PAGE

Anti-TNFRSF10A / TRAILR1 / CD261 Antibody (mapatumumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



SEC-HPLC

The purity of Anti-TNFRSF10A / TRAILR1 / CD261 Antibody (mapatumumab) is 98.89%, determined by SEC-HPLC.