

# Human LDHA/PIG19 (His tag) recombinant protein



Catalog Number: 518143

## General Information

### Gene Name Synonym

L-Lactate Dehydrogenase A Chain, LDH-A, Cell Proliferation-Inducing Gene 19 Protein, LDH Muscle Subunit, LDH-M, Renal Carcinoma Antigen NY-REN-59, LDHA

### Protein Construction

A DNA sequence encoding the Human LDHA/PIG19 Met1-Phe332 is expressed with a 6His tag at the N-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 100mM NaCl, 50% Glycerol, pH 8.0.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human LDHA/PIG19 (His tag) recombinant protein