

Human LDHA/PIG19 (His tag) recombinant protein



Catalog Number: 518143

General Information

Gene Name Synonym

L-Lactate Dehydrogenase A Chain, LDH-A, Cell Proliferation-Inducing Gene 19 Protein, LDH Muscle Subunit, LDH-M, Renal Carcinoma Antigen NY-REN-59, LDHA

Protein Construction

A DNA sequence encoding the Human LDHA/PIG19 Met1-Phe332 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM TrisHCl, 100mM NaCl, 50% Glycerol, pH 8.0.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human LDHA/PIG19 (His tag) recombinant protein