

VPS37A antibody



Catalog Number: 107543

Product name

VPS37A antibody

Immunogen

[Human VPS37A Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

VPS37A Mouse Monoclonal antibody. Positive WB detected in SW 1990 cells, MCF-7 cells. Observed molecular weight by Western-blot: 44 kDa

Preparation

This antibody was obtained by immunization of VPS37A recombinant protein (Accession Number: XM_017013026). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:1000-1:10000

Validations



SW 1990 cells were subjected to SDS PAGE followed by western blot with Catalog No:107543(VPS37A Antibody) at dilution of 1:2000