

Anti-Estrogen Receptor alpha (1G4) Mouse antibody



Catalog Number: 168543

Product name

Anti-Estrogen Receptor alpha (1G4) Mouse antibody

ICC: 1/200 - 1/1000

ELISA: 1/10000

Validations

Specificity

Human

Antibody description

Estrogen Receptor alpha (1G4) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human ESR1 expressed in E. Coli.

Formulation

Ascitic fluid containing 0.03% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

WB, IHC-P, ICC/IF


Dilutions

WB: 1/500 - 1/2000


IHC: 1/200 - 1/1000



Immunohistochemical analysis of paraffin-embedded mammary cancer tissues using ESR1 mouse mAb with DAB staining.



Immunofluorescence analysis of HeLa cells using ESR1 mouse mAb (green). Red



Western blot analysis of extracts from SHSY-5Y and Ramos cell lysates using Estrogen Receptor α mouse mAb (1:500 diluted). Predicted band size: 66 kDa. Observed band size: 66 kDa.